

Contents

1	Introduction to the system and its purpose (outline briefly what it does).	2
2	How to install the program/database on your system.	3
2.1	How to install	4
3	What to expect as you enter the system (password, menu etc)	9
4	How to use the system:	10
4.1	The menus	10
4.1.1	Calculators	10
4.1.2	File	11
4.2	Explain what each form does	11
4.3	Calculators	12
4.3.1	Body Surface Area (BSA)	12
4.3.2	Validation	13
4.3.3	Renal Function	14
4.3.4	Calvert	15
4.4	Regime Viewer & Editor	16
4.4.1	How to add a new Regime	16
4.5	Regime	18
4.6	Comment on any validation the user can expect to happen . .	18
4.7	Comment on any on-screen help/menus etc	18
4.8	Give samples of the error messages generated (in-situ!)	18
4.9	Indicate what to do if something goes wrong	18
5	Explain how to back-up the system (program and data files).	19

Chapter 1

Introduction to the system and its purpose (outline briefly what it does).

This is a software aid for pharmacists who need to prepare and dispense Chemotherapy Regimes for patients being treated for neoplasms.

This software uses an internal user modifiable record file to store the regime data. The software is also records and stores patient data, using it to calculate the required doses on the appropriate days.

Chapter 2

How to install the program/database on your system.

On a company network it is recommended that the installation is left to IT support to ensure that the files are placed in central location that is accessible to all who require it.

This also ensures that it can be made part of the system wide backup procedure.

The instructions below cover the installation method for a single user on Windows 7

2.1 How to install

1. Get the Zipped Regime file, right click on it and choose **Extract All**.




Figure 2.1

2. When this window appears click **Browse**.




Figure 2.2

3. To place the file in your Start Menu, Click on your User name then AppData ⇒ Roaming ⇒ Microsoft ⇒ Windows ⇒ Start Menu ⇒ Programs, now click Make New Folder and name it Regime. Press Enter and click Ok.




Figure 2.3

4. The Filepath should now look like this C:\Users\A.2\AppData\Roaming\Microsoft\Windows\Start Menu\Programs\Regime




Figure 2.4

5. Click Extract. Providing the Show extracted files when complete tickbox was selected a file explorer window will open with the selected folder.




Figure 2.5

You may need to log out for the file to be added to the start menu, but it is possible to search for it.




Figure 2.6

Once it is visible it is possible to pin it to the start menu or the Taskbar by using the options from the right click menu.




Figure 2.7

Chapter 3

What to expect as you enter the system (password, menu etc)

When the program opens you will see the interface shown below.




Figure 3.1

Clicking on **Calculators** will call up a self contained calculator interface, useful for quick calculations.

Clicking on **Regime** will open the main regime dosage system.

Under the File Menu is the Regime Editor, this is used to view and edit the data on Regimes.

Chapter 4


How to use the system:

4.1 The menus

4.1.1 Calculators

This menu calls the Calculators form, taking you directly to the appropriate tab. These options can also be triggered by using the keyboard shortcuts, these work on all forms that include this menu. The Shortcuts are:

- Alt-B BSA
- Alt-C Cockcroft-Gault
- Alt-A Calvert



4.1.2 File

The Regime Editor is used to review and if necessary edit the regime record file.

4.2 Explain what each form does

Start

Gives you access to the Calculators, Regime and Regime Viewer & Editor.

Calculators

Self-contained calculators for BSA, Renal Function and Calvert.

Regime

Displays a list of patients, includes search and filter by regime functionality.

PatientData

This form is used to view and update patient data and includes integrated dose calculators.

RegimeViewer

Displays a list of Regimes.

RegimeEditor

This form is used to view and update Regimes.

SPC Viewer

Specification of Product Characteristics viewer, this shows the SPC page for the current drug as shown on the website

4.3 Calculators

This form allows you to quickly calculate values required for determining suitable regime doses. The included calculators are for BSA (Dubois and Dubois method), the Creatinine Clearance level (Cockcroft-Gault method) and the Calvert Formula.

You can access this form via the **Calculator** Button on the start page, or from the **Calculator Menu**.

4.3.1 Body Surface Area (BSA)


This system uses the Dubois & Dubois method with the formula

$$BSA = 0.20247 \times Height[m]^{0.725} \times Weight[Kg]^{0.425}$$

The BSA calculator requires the patients height in meters and the weight in Kilograms; this produces a value in meters squared.

BSA limit

Due to the increased risk, current guidelines recommend that BSA values are capped at 2.0. If the BSA is more than 2.0 you will see this warning:




You should discuss the case with the Consultant to determine if a higher dose is suitable in the circumstances.


4.3.2 Validation

Both inputs accept decimal values greater than zero, any inputs that don't fulfill this criteria (zero, negative values or text) will result in an error message and the relevant input field will have a red background.

Once you have a valid answer you can click on the **Use result for Calvert** to switch directly to the Calvert Calculator with the eCCR textbox already filled. Note: this button is not visible until a valid result has been calculated for the eGRF



(a) Valid Input



(b) Invalid Input

Figure 4.1: BSA Calculator in use

4.3.3 Renal Function

The Cockcroft-Gault formula is used to produce the estimated Glomerular Filtration Rate (eGFR).

$$eGFR[\text{ml/min}] = \frac{(140 - \text{age}[\text{years}]) \times \text{Mass}[\text{Kg}] \times k(\text{GenderConstant})}{\text{SerumCreatinine}[\mu\text{mol / L}]}$$

$$k_{\text{Male}} = 1.32$$

$$k_{\text{Female}} = 1.04$$

The year is to be given as an integer, Mass and Serum Creatinine are decimals and the gender constant is set by selecting the tick box for male patients.


(a) Valid Input

(b) Invalid Input

Figure 4.2: eGFR Calculator in use


4.3.4 Calvert

The Calvert Calculator requires that both the eGFR and the Area Under the Curve (AUC) are integers with a value greater than or equal to 1.



The screenshot shows a window titled "Calculators" with three tabs: "BSA", "Cockcroft-Gault", and "Calvert". The "Calvert" tab is selected. It contains three input fields: "Target AUC" with the value "5", "eGFR [ml/min]" with the value "34", and "Dose [mg]" with the value "295". Below the input fields are two buttons: "Calculate Dose" and "Get CCR".

(a) Valid Input



The screenshot shows the "Calculators" window with the "Calvert" tab selected. The "Target AUC" field contains the value "5", which is highlighted in red. Below it, the "eGFR [ml/min]" field contains the value "34". The "Dose [mg]" field is empty. An "Input Error" dialog box is open in the foreground, displaying a red "X" icon and the text: "Out of Range Error" and "The Given value is either too big or too small for this field". An "OK" button is at the bottom right of the dialog box.

(b) Invalid Input

Figure 4.3: BSA Calculator in use

4.4 Regime Viewer & Editor

Accessed via the File Menu on the start page or the keyboard shortcut **Alt-R**, this interface is used to view and if required edit the data for regimes. This shows the form when it is loaded for the first time if no initial data is provided.





Figure 4.4

4.4.1 How to add a new Regime

1. Click on **New Regime** (bottom right hand corner).
2. This opens the **RegimeEditor** form as shown below.



The image shows a software window titled "RegimeEditor". At the top, there is a header bar with a green icon and the title. Below the header, there are three input fields: "RegimeName" (a text box), "#Cycles" (a spinner box with the value 0), and "Days/Cycle" (a spinner box with the value 0). Below these fields is a tabbed interface with two tabs: "tabPage1" and "Dosage". The "Dosage" tab is currently selected. Below the tabs, there are three text boxes: "Extravasation" and "Comment" (which is larger and occupies more space). At the bottom right of the window, there is a "Save" button.

Figure 4.5

Using a reputable source of Regime data copy and paste in the data as follows:

RegimeName Enter the name of the Regime, this usually takes the form of an acronym derived from the Active Ingredients.

#Cycles The standard number of cycles of a regime that a patient will undergo.

Figure 4.6

Days/Cycle The number of days per cycle.

First textbox on the first tab A quick summary of the Regime, including its purpose.

Extravasation Discusses the risk of patients undergoing extravasation while on this regime.

Comment This is a section for additional data including tests to be performed and modifications that may be required.

4.5 Regime

The main functionality of this software is accessed via the **Regime** button. Once you have clicked on this button you will see the form shown below. The list on the left displays the names of all the available regimes, the list on the right shows all the current patients. The textbox under the patient list is a search box, enter any part of the patients name to filter the list of patients. While the system should start searching automatically, there is also a search button.

4.6 Comment on any validation the user can expect to happen

4.7 Comment on any on-screen help/menus etc

4.8 Give samples of the error messages generated (in-situ!)

4.9 Indicate what to do if something goes wrong

Chapter 5

Explain how to back-up the system (program and data files).